

# CytoFlexLX

CytoflexLX A and B have **almost** the same configurations. Only one detector difference off violet laser. VSSC vs. BV786

	<b>Laser (NM)</b>	<b>Bandpass</b>	<b>Channel Name</b>	<b>Additional Optimal Colors</b>
LASER	355	405/30	<b>BUV395</b>	
		450/45	<b>DAPI</b>	ZombieUV, L/Dfix Blue, DyLight 350
		740/35	<b>BUV737</b>	
LASER	405	405/10	<b>VSSC (CytoA only)</b>	Violet Side Scatter - used for small particle detection
		450/45	<b>PacificBlue</b>	BV421, V450, Alexa405, CF405M, DyLight405, eFluor450, SuperBright436, VioBlue, Zombie Violet, L/DfixViolet
		525/40	<b>BV510</b>	L/Dfix Aqua, V500, AF430, DyLight405LS, VioGreen, Zombie Aqua
		610/20	<b>BV605</b>	Qdot605, DyLight405LS, eFluor605NC, eVolve605, SuperBright600
		712/25	<b>BV711</b>	Qdot705, SuperBright702
		763/43	<b>BV786 (CytoB only)</b>	Qdot800
LASER	488	525/40	<b>FITC</b>	GFP, Alexa488, CFSE, BB515, DyLight488, Zombie Green, L/Dfix Green
		610/20	<b>B610-ECD</b>	<i>all current available fluors for this channel are best excited by the 561 laser</i>
		690/50	<b>PerCP</b>	PerCP-Cy5.5, BB700, PerCP-eFluor710, PerCP-Vio700
LASER	561	584/42	<b>PE</b>	dsRed, RFP
		610/20	<b>PE-TexasRed</b>	mCherry, PE-Dazzle, PE-CF594, PE-Vio615, PE-eFluor610, DyLight594, Zombie Red, L/Dfix Red
		675/30	<b>PE-Cy5</b>	mPlum
		710/50	<b>PE-Cy5.5</b>	
		763/43	<b>PE-Cy7</b>	PE-Vio770
LASER	638	660/10	<b>APC</b>	AF647, Dylight633, L/Dfix Far Red
		712/25	<b>AL680</b>	Alexa700
		763/43	<b>APC-Cy7</b>	APC-Fire750, APC-Alexa750, APC-eFluor780, APC-Vio770, Zombie NIR, L/Dfix NIR

\* **L/Dfix** = Live/Dead Fixable

Please ask Facility Staff if you have questions or concerns about particular fluorochromes or fluorescent proteins