


Yale CryoEM Resource

Tecnai T12 or Talos L120C Access Request ^a

User and billing information		
1	User name	Jane Doe
2	E-mail	Jane.doe@yale.edu
3	PI name(s)	Joe Bloggs
4	E-mail(s)	Joe.bloggs@yale.edu
5	COA code	GR022222.CC0222.PG00032.PJ000001.netid
Scheduling preference		
6	Instrument (YSM Tecnai T12 or SH Talos L120C) ^b	SH Talos L120C
7	Favorite start/end date and time	7/10, 10 am-12 pm
8	Alternative dates and time	7/15, 10 am-12 pm
Sample information		
9	Experiment type (A. Negative Stain EM; B. Cryo sample screening; C. TEM imaging; D. Tomography)	A
10	Grid arrival date	7-2-2025
11	Number of grids	8
12	 Grid Box Labels and Grid position (Make your labels LEGIBLY)	Box-1, Position 1-3 Box-2, Position 2, 3

^a Refer to the **YCR Policies** document for allowed samples and sign-up rules for each instrument; **This form is only required for requesting Attended sessions in which YCR will provide partial or full service.**

^b YSM Tecnai T12 is located in SHM IE room 26E; SH Talos L120C is located in Bass Center Room 231.

Preliminary data (optional, but are beneficial for the experiments)

- I. **Tomography:** Provide a brief description of your sample and target of interest, and indicate if you plan to use subtomogram averaging.
- II. **SPA:** It is recommended to include images demonstrating sample purity or quality, such as an SDS-PAGE gel image or previous negative-stain EM micrograph.